
INSTRUCTIONS FOR USE







EN	Instructions for use	3
CS	Návod k použití	7
DA	Brugervejledning	11
DE	Gebrauchsanweisung	15
ES	Instrucciones de uso	19
FI	Käyttöohje	23
FR	Notice d'utilisation	27
IT	Istruzioni d'uso	31
NL	Gebruiksaanwijzing	35
NO	Bruksanvisning	39
PL	Instrukcja użytkowania	43
PT	Instruções de utilização	47
RO	Instrucțiunile de utilizare	51
RU	Инструкция по применению	55
SV	Bruksanvisning	59

STRUCTURE OF DIPSTICK

TEST PROCEDURE AND RESULTS

STRUCTURE OF CARTRIDGE

Instructions for use

Numbers 1 - 13 refer to illustrations on inner cover.

Structure of dipstick

1 Dip area 2 Result area 3 Test line 4 Control line

Structure of cartridge

5 QR code

Intended use

Actim® Partus 1ngeni test is a quantitative immunochromatographic dipstick test for measuring the concentration level of phosphorylated IGFBP-1 (insulin-like growth factor binding protein-1) in cervical secretions during pregnancy. The test is read using the Actim 1ngeni instrument. The test is intended for professional use to help predict the risk of preterm or imminent delivery when fetal membranes are intact. A low risk test result (pHIGFBP-1 concentration less than 10 µg/l) is a clear indication that the patient will not deliver within 7-14 days.

Kit components

Actim Partus 1ngeni kit (31931RETAC) contains 10 test packs with instructions for use. The components of each Actim Partus 1ngeni test pack (31921RETAC) are:

- One sterile polyester swab for specimen collection.
- One tube of Specimen Extraction Solution (0.5 ml). This phosphate-buffered solution contains bovine serum albumin (BSA), protease inhibitors and preservatives.
- One dipstick in a sealed aluminum foil pouch with desiccant.
- One cartridge for reading the result with the Actim 1ngeni instrument. The QR code on the cartridge includes information about the test (Actim Partus 1ngeni), the dipstick lot, expiry date, and lot specific calibration data.

Storage

Store the test kit at +2... +25° C. Stored unopened, each component can be used until the expiry date marked on the component. The kit can also be stored for 2 months at +2... +30° C. Use the dipsticks shortly after their removal from the aluminum foil pouch.

Specimen collection and sample handling

The specimen consists of cervical secretion that is extracted into the Specimen Extraction Solution provided. A sample is obtained using a sterile polyester swab (provided in the test pack).

1 The sample should be collected prior to performing a digital examination and/or transvaginal ultrasound. 1 Take care not to touch anything with the swab before taking the sample.

1. A cervical secretion sample is obtained using a sterile polyester swab (provided in the kit) from the cervical os during a sterile speculum examination.
2. The swab should be left in the cervical os for 10-15 seconds to allow it to absorb the secretion specimen 6.
3. Open the Specimen Extraction Solution tube and put it in a vertical position. The specimen is extracted immediately from the swab by swirling the swab vigorously in the extraction solution for 10-15 seconds 7. Press the swab against the wall of the Specimen Extraction Solution tube to remove any remaining liquid from the swab. Discard the swab.

1 Specimens should be tested as soon as possible after extraction but, in any case, no more than 4 hours after the specimen collection and extraction. If a specimen cannot be tested within this time it should be frozen. After thawing, the specimens should be mixed and tested as described below.

Test procedure – Patient Test

Test reaction occurs inside the instrument and results are read automatically at 5 minutes. Insertion of the cartridge into the Actim 1ngeni automatically initiates the Actim Partus 1ngeni test procedure.

1 If the test kit is stored refrigerated, allow the test pack to reach room temperature.

Prepare the instrument:

1. Switch the Actim 1ngeni on (if power off). After the automatic self test, log onto the instrument by entering your user ID either manually or with a barcode reader.
2. Select Test icon 8 and select Patient Test 9.
3. Enter patient ID either manually or with a barcode reader.

Activate the test and start the procedure:

1. When the Actim 1ngeni screen shows “Activate and insert test”, open the foil pouch containing the dipstick.

1 Do not touch the yellow dip area at the lower part of the dipstick. Identifying marks may be written on the upper purple part of the dipstick. The dipstick must be used shortly after its removal from the foil pouch.

2. Activate the test by placing the yellow dip area of the test into the extracted sample **10** and hold it there until you see the liquid front enter the result area **11**.
 3. Immediately remove the dipstick from the solution, place it in the cartridge (reading area forward and coloured side facing up) and insert the cartridge inside the Actim 1ngeni (QR code forward and facing up) until the cartridge clicks into place **12**. Follow and respond to any information on the screen.
- !** Use only the cartridge included in the test pack. **!** Ensure that the dipstick is firmly inside the cartridge. **!** Ensure that the timer appears on the screen.

Proceed to section "Results".

Results

The Actim 1ngeni automatically reports and saves the result. Result is reported as a **low risk** result when the phIGFBP-1 concentration is less than 10 µg/L and as an **elevated risk** result when the phIGFBP-1 concentration is equal or more than 10 µg/L.

- !** Do not visually interpret the result under any circumstances after removing the test from the instrument because the test is optimised for the Actim 1ngeni.
- When the test result is visible on the screen, remove and discard the dipstick and the cartridge. The result can be either printed or exported into a USB stick or stored in the Actim 1ngeni for later use. Refer to the instrument manual for details.
 - Log off when leaving the instrument.

Test result with additional test-related information are shown on the screen **13**. If problems have occurred during the test process, an error is reported and no result is shown.

Quality control and Calibration

Test calibration and quality control

Each Actim Partus 1ngeni test kit lot has been calibrated and quality controlled during manufacturing to give reproducible results with the Actim 1ngeni instrument. QR code on the cartridge contains the lot specific calibration data. Depending on local hospital practices and the local administrator for the Actim 1ngeni instrument, quality control for the Actim Partus 1ngeni test system using liquid controls may be advised. Refer to Actim Partus Controls (31900ETAC) for details. There are no daily requirements for controlling the instrument performance using liquid controls.

Actim 1ngeni instrument

The Actim 1ngeni instrument has various levels of quality control procedures available. Refer to Actim 1ngeni instrument manual for details.

Limitations of the test

- The test is intended for *in vitro* diagnostic use only.
- Before performing the test ensure that the fetal membranes are intact (for example with Actim PROM test), because with ruptured fetal membranes Actim Partus test will also give a positive result.
- Patients with moderate or heavy vaginal bleeding should not be tested. It is recommended that a sample be taken when bleeding has stopped, and the extract is essentially blood-free.
- The test result indicates the risk at the time of sampling and changes in the patient's condition may later affect the final outcome of the pregnancy.

Notes

- The test requires approximately 150 µl of extracted sample to ensure proper performance of the test.
- Use care when placing the dipstick in the sample tube. The upper part of the dipstick must stay dry.
- Do not use a dipstick that has become wet before use because moisture damages the dipstick.
- Do not use a dipstick if you notice a blue coloring in the result area before testing.
- Do not use the dipstick if its aluminum foil pouch or the seals of the pouch are not intact.
- Use only the swab provided with the Actim Partus kit.
- Improper sampling may lead to false result.
- When dipping, be careful to hold the dipstick in position (with the dip area in the sample extract) until the sample liquid front reaches the result area.
- Do not insert anything other than the cartridges provided with Actim 1ngeni tests into the instrument.
- Proper timing of reading the dipstick is essential. A delay in inserting the dipstick into the instrument can have an effect on the test results. In these cases the test should be repeated.
- Ensure that the dipstick is firmly inside the cartridge and does not slip out when placed into the reader.

- If the procedural control line is not detected by the Actim 1ngeni instrument, the result will be declared indeterminate and the test should be repeated.
- The result must not be visually inspected after analysing it with the instrument.
- Do not re-analyse a test with the instrument.
- For additional information, please refer to the Actim 1ngeni instrument manual.
- As with all diagnostic tests, results must be interpreted in the light of other clinical findings.
- All biological specimens and materials must be treated as potentially hazardous, and disposed of in accordance with local authority guidelines.

Principle of the test

Decidual cells synthesize the phosphorylated forms of IGFBP-1 (pHIGFBP-1) while amniotic fluid contains substantial quantities of non- and less phosphorylated forms of IGFBP-1. When delivery is approaching, fetal membranes begin to detach from the decidua parietalis, and small amounts of pHIGFBP-1 begin to leak into cervical secretions. In Actim Partus 1ngeni test a cervical specimen sample is taken with a sterile polyester swab during sterile speculum examination and extracted into the Specimen Extraction Solution. The quantity of pHIGFBP-1 in the solution is detected using a dipstick and Actim 1ngeni instrument.

The test is based on immunochromatography. It uses two monoclonal antibodies to human IGFBP-1. One is bound to blue latex particles [the detecting label]. The other is immobilized on a carrier membrane to catch the complex of antigen and latex-labeled antibody and indicate a positive result. When the dip area of the dipstick is placed in an extracted sample, the dipstick absorbs liquid, which starts to flow up the dipstick. If the sample contains pHIGFBP-1 it binds to the antibody labeled with latex particles. The particles are carried by the liquid flow and, if pHIGFBP-1 is bound to them, they bind to the catching antibody. A blue line (test line) will appear in the result area if the concentration of pHIGFBP-1 in the sample exceeds the detection limit of the test. A second blue line (control line) confirms correct performance of the test. The lines are read, and intensity of line color is quantitated with Actim 1ngeni instrument.

Performance characteristics

Clinical performance study 1

The clinical performance of the Actim Partus 1ngeni test was evaluated in a retrospective study using specimens obtained from pregnant women with

symptoms of preterm delivery (n=58; GA 22-42 weeks) and controls (n=62; GA 22-42 weeks).

Assessment of quantitative test results was made using delivery within 7 days as an end point. Comparison was also made to the cervical length. The concentration of pHIGFBP-1 was higher in those women, who delivered within 7 days after sample collection. Also, cervical canal was shorter in those patients. Statistical analysis is shown in Table. 1 on the inner back cover. Correlation between clinical outcome and concentration of pHIGFBP-1, cervical length and the combination of these are statistically significant (p-value < 0.05) when delivery occurs within 7 days (Table. 1).

Clinical performance study 2

The clinical performance of the quantitative Actim Partus 1ngeni was assessed for pregnant women with signs and symptoms of preterm labor. This retrospective study included 691 samples obtained from symptomatic women (GA 22-35) with singleton pregnancies. From this group 360 woman delivered spontaneously.

The clinical performance and statistical analyses showed, that the increasing concentration of pHIGFBP-1 correlated with increased relative risk of delivery within 7 and 14 days from sample collection. Similarly, an increasing concentration of pHIGFBP-1 correlated with increased risk of spontaneous labor and preterm birth before 30, 34 and 37 weeks of gestation. These findings and their statistical significance are summarized in Table 2, 3, 4 and 5.

Analytical sensitivity

Analytical sensitivity (detectability) was tested with the quantitative Actim Partus 1ngeni test. Negative samples as well as low level samples were used for testing. The limit of blank, LoB, of the Actim Partus 1ngeni test was found to be approximately 0.7 µg/L. The detection limit, LoD, of the quantitative Partus test was found to be approximately 1.5 µg/L. The limit of quantification (LoQ) was through precision study found to be approximately 3.5 µg/L.

Analytical specificity

Analytical specificity (cross-reactivity) was tested with human IGFBP proteins. The IGF binding proteins -2, -4, -5 and -6 were tested at a concentration level of 5000 µg/l and IGF binding protein -3 at a concentration level of 50,000 µg/l. No interference of the IGF binding proteins was observed with the performance of Actim Partus 1ngeni test. The test is specific to human IGFBP-1.

Repeatability

A repeatability study was performed. Three panels of specimens were tested on

twenty days with two runs per day and two replicates per run (80 replicates per sample). Sample concentrations and results are shown in Table 6.

Reproducibility

A Reproducibility study was performed. Six panels of specimens were tested on five days with one run per day and five replicates per run using three instruments (75 replicates per sample). The repeatability, within-laboratory and reproducibility precision coefficients of variation, CV% are shown in Table 7.

Interference substances

The following substances were tested with the Actim Partus 1ngeni test and were found not to affect the performance of the test.

INTERFERING SUBSTANCE	CONCENTRATION TESTED
Surgical lubricant	50%
Betadine	50%
Econazole cream	250 mg/mL
Hexachlorophene cream	250 mg/mL
Progesterone capsule	250 mg/mL
Ultrasound transmission gel	10%
Personal Lubricant	10%
Baby powder	50%
Baby oil	50%
Feminine Deodorant Film	50%
Vaginal Gel	10%
Shower and bath products (various)	1%
Bilirubin	500 µg/mL

- Pregnancy urine was tested with the Actim Partus 1ngeni test. No interference of pregnancy urine was observed with the performance of the test.
- Amniotic fluid was tested with the test. It was observed, that amniotic fluid interferes with the performance of the test.
- Whole blood was tested with the test. It was observed, that whole blood interferes with the performance of the test.

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Table 1.

CORRELATION BETWEEN CLINICAL OUTCOME AND CLINICAL DATA WITHIN 7 DAYS.

Outcome	phIGFBP-1 median ($\mu\text{g/l}$)	Cervical length median (mm)	p-value
≤ 7 days	15.0	15.0	-0,421*
> 7 days	4.0	30.0	0.076
p-value	0.011*	0.001*	*p ≤ 0.05

Table 3.

RELATIVE RISK OF PRETERM DELIVERY ACCORDING TO QUANTITATIVE ACTIM PARTUS 1NGENI CATEGORIES.

phIGFBP-1 category (ng/ml)	Delivery ≤ 7 days	Delivery ≤ 14 days	Delivery < 30 wks	Delivery < 34 wks	Delivery < 37 wks
≤ 2	-	-	-	-	-
2.1-9.9	5.3*	1.5	0.5	1.4	1.6*
10.0-49.9	5.8*	2.1	5.8*	3.3**	2.4***
50.0-249.9	7.7*	4.2*	11.5**	4.1**	2.4***
≥ 250	13.5**	11.1***	20.3***	8.7***	3.9***

Relative risk compared to phIGFBP-1 ≤ 2 ng/ml.z-statistics significance: *P <0.05 , **P <0.01 , ***P < 0.001

Table 2.

PRETERM DELIVERY RATES IN SYMPTOMATIC WOMEN ACCORDING TO QUANTITATIVE ACTIM PARTUS 1NGENI CATEGORIES.

phIGFBP-1 category (ng/ml)	Delivery ≤ 7 days	Delivery ≤ 14 days	Delivery < 30 wks	Delivery < 34 wks	Delivery < 37 wks
≤ 2	0.7%	2.6%	0.7%	3.3%	15.7%
2.1-9.9	3.8%	3.8%	0%	4.5%	24.8%
10.0-49.9	4.1%	5.4%	4.1%	10.8%	37.8%
50.0-249.9	5.4%	10.8%	8.1%	13.5%	37.8%
≥ 250	9.5%	28.6%	14.3%	28.6%	61.9%

Table 4.

SPONTANEOUS PRETERM DELIVERY RATES IN SYMPTOMATIC WOMEN ACCORDING TO QUANTITATIVE ACTIM PARTUS 1NGENI CATEGORIES.

phIGFBP-1 category (ng/ml)	Delivery ≤ 7 days	Delivery ≤ 14 days	Delivery < 30 wks	Delivery < 34 wks	Delivery < 37 wks
≤ 2	0.5%	3.4%	1.0%	3.9%	19.5%
2.1-9.9	2.9%	2.9%	0%	5.7%	28.6%
10.0-49.9	4.9%	7.3%	4.9%	9.8%	39.0%
50.0-249.9	7.7%	15.4%	11.5%	19.2%	42.3%
≥ 250	11.1%	27.8%	16.7%	27.8%	61.1%

Table 5.

RELATIVE RISK OF SPONTANEOUS PRETERM DELIVERY ACCORDING TO QUANTITATIVE ACTIM PARTUS 1NGENI CATEGORIES

phIGFPB-1 category (ng/ml)	Delivery ≤ 7 days	Delivery ≤ 14 days	Delivery < 30 wks	Delivery < 34 wks	Delivery < 37 wks
≤2	-	-	-	-	-
2.1-9.9	5.9	0.8	0.6	1.5	1.5
10.0-49.9	10.0	2.1	5.0	2.5	2.0**
50.0-249.9	15.8*	4.5*	11.8**	4.9**	2.4**
≥ 250	22.8**	8.1***	17.1**	7.1***	3.1***

Relative risk compared to phIGFPB-1 <2 ng/ml.

z-statistics significance: *P<0.05, **P<0.01, ***P < 0.001

Table 6.

RESULTS FROM THE REPEATABILITY STUDY.

Sample	Mean value, µg/l	Repeatability		Within-Laboratory Precision (overall)	
		SD, µg/l	%CV	SD, µg/l	%CV
1	3.5	0.46	13.1	0.69	19.7
2	14	1.50	10.7	2.33	16.6
3	43	4.87	11.3	6.01	14.0

Table 7.

RESULTS FROM THE REPRODUCIBILITY STUDY.

Sample	Mean value, µg/l	N	Repeatability		Between-Run*		Between-Day		Between-Instrument		Reproducibility	
			SD, µg/l	%CV	SD, µg/l	%CV	SD, µg/l	%CV	SD, µg/l	%CV	SD, µg/l	%CV
1	3.5	75	0.61	17.4	NA	NA	0.85	24.3	0.36	10.3	1.10	31.4
2	15	75	2.16	14.4	NA	NA	2.25	15.0	1.97	13.1	3.69	24.6
3	40	75	7.21	18.0	NA	NA	4.41	11.0	4.16	10.4	9.42	23.6
4	157	75	22.29	14.2	NA	NA	25.91	16.5	13.00	8.3	36.57	23.3
5	341	75	57.52	16.9	NA	NA	78.48	23.0	13.59	4.0	98.25	28.8
6	778	75	139.0	17.9	NA	NA	262.5	33.7	NA	NA**	297.0	38.2

* One run per day performed

** A negative variance

Literature • Literatura • Litteratur • Literatur • Bibliografia • Kirjallisuus • Bibliographie • Bibliografia • Literatuur • Literatur • Literatura • Literature • Literatură de specialitate • Литература • Litteratur

Akercan F et al.

Value of cervical phosphorylated insulinlike growth factor binding protein-1 in the prediction of preterm labor. *J Reprod Med* [2004] 49: 368-372.

Altinkaya O et al.

Cervical phosphorylated insulin-like growth factor binding protein-1 in prediction of preterm delivery. *Arch Gynecol Obstet.* [2009] 279:279-283.

Azlin MI et al.

Role of pIIGFBP-1 and ultrasound cervical length in predicting pre-term labour. *Journal of Obstetrics and Gynaecology* [2010] 30:456-460.

Brik Spinelli M et al.

Phosphorylated insulin-like growth factor binding protein-1 and cervical measurement in women with threatening preterm birth. *Acta Obstet Gynecol Scand* [2010] 89:268-74.

Eroglu D et al.

Prediction of preterm delivery among women with threatened preterm labor. *Gynecol Obstet Invest* [2007] 64:109-116.

Kekki M et al.

Insulin-like growth factor-binding protein-1 in cervical secretion as a predictor of preterm delivery. *Acta Obstet Gynecol Scand* [2001] 80:546-551.

Lembet A et al.

New rapid bed-side test to predict preterm delivery: phosphorylated insulin-like growth factor binding protein-1 in cervical secretions. *Acta Obstet Gynecol Scand* [2002] 81:706-712.

Rutanen EM

Insulin-like growth factors in obstetrics. *Opin Obstet Gynecol* [2000] 12:163-168.

Tanir HM , Sener T, Yildiz Z.

Cervical phosphorylated insulin-like growth factor binding protein-1 for the prediction of preterm delivery in symptomatic cases with intact membranes. *J Obstet Gynaecol Res* [2009] 1:66-72.

Ting HS et al.

Comparison of bedside test kits for prediction of preterm delivery: phosphorylated insulin-like growth factor binding protein-1 (pIGFBP-1) test and fetal fibronectin test. *Ann Acad Med Singapore* [2007] 36:399-402.

Actim[®] Partus 1ngeni

Actim Oy
Klovinpellontie 3, FI-02180 Espoo, Finland
Tel. +358 9 547 680
www.actimtest.com
actim@actimtest.com

© Actim Oy. All rights reserved.
Actim is a trademark of Actim Oy.

AORACE31931-2 2020-05-14